

| AUDIOMETERS | | | | | | | | | | | | | | | |
|--|----------------------|------------------------|------------------------------|-----------|--------------------|----------------------|-----------------|---------------------------|---------------------|--|---------------------|----|------------------------|----------|----------------------|
| Manufacturer | Frequency Range (Hz) | Attenuation Range (db) | Test Program/Control Options | | | Test Results Options | | | Time Per Test (min) | # of Persons Tested Concurrently From One Unit | Is a Booth Required | | Audiometer Calibration | | # of Validity Checks |
| | | | Manual | Automatic | Computer Automated | Printout | Database Stored | Automatic Data Processing | | | Yes | No | Internal | External | |
| Beltone Electronics Corp. 773/583-3600 | 125-8000 | -10-110 | • | | | | | | 1 | | • | | • | | |
| Maico Diagnostics 888/941-4201 | 500-8000 | -10-100 HL | • | • | • | • | • | • | 6 | 1-8 | | • | | • | |
| Micro Audiometrics Corp. 828/644-0771 | 125-8000 | -10-100 | • | • | • | • | • | • | 3-5 | 1 | | • | | • | |
| Occupational Health Dynamics 205/980-0180, 888/464-3872 | 500-8000 | 0-100 | • | • | • | • | • | • | 3-5 | 16 | | • | | • | |
| Tremetrics Occ. Health Group 800/825-0121 | 500-8000 | -10-100 | • | • | • | • | • | • | 3-6 | 8 | | • | • | • | |